


Microscopy Sources

Product Catalog



LIGHT CONVERSION is a global leader in ultrafast technology, designing and manufacturing:

- > Femtosecond Lasers,
- > Wavelength-Tunable Sources,
- > OPCPA Systems,
- > Microscopy Sources,
- > Spectroscopy Systems.

The comprehensive portfolio represents the best-in-class lasers tailored for industry, science, and medicine.

About Us

Founded in 1994, LIGHT CONVERSION has evolved into a leading company in ultrafast laser technology with over 9000 systems installed worldwide and 600 employees, 15% of whom focus on R&D. The company's lasers are used in all of the top 50 universities worldwide, highlighting its commitment to state-of-the-art research, while also ensuring the reliability and performance in 24/7 industrial applications. With international offices in the US, China, and Korea, along with a global representative network, the company ensures worldwide sales and service.

Contents

5 Microscopy Sources

CRONUS | 2P

- 6 Three-Channel Wavelength-Tunable Femtosecond Laser

CRONUS | 3P

- 8 Laser Source for Advanced Nonlinear Microscopy

11 Wavelength-Tunable Sources

I-OPA

- 12 Industrial-Grade Optical Parametric Amplifier

ORPHEUS | F

- 14 Broad-Bandwidth Hybrid Optical Parametric Amplifier

ORPHEUS | TWINS

- 16 Dual Optical Parametric Amplifier

19 Femtosecond Lasers

FLINT

- 20 High-Repetition-Rate Lasers

PHAROS

- 22 Modular-Design Femtosecond Lasers for Industry and Science

CARBIDE

- 26 Unibody-Design Femtosecond Lasers for Industry and Science

31 Nonlinear Microscopy Applications

36 Global Representative Network





Microscopy Sources

The CRONUS femtosecond lasers cover applications in functional neuroimaging, optogenetics, and deep imaging, using medium-repetition-rate three-photon (3P) excitation and fast high-repetition-rate two-photon (2P) imaging, as well as widefield and holographic excitation.

CRONUS | 2P

Three-channel laser with a high repetition rate for simultaneous 2P excitation of multiple fluorescent probes, calcium indicators, opsins or CARS and SRS.

CRONUS | 3P

Turn-key laser source with μJ -level pulses, covering the biological transparency windows at 1300 and 1700 nm for 3P microscopy and 1030 nm for optogenetic stimulation.

Optimized for advanced multiphoton microscopy

Plug-and-play functionality with automated wavelength and dispersion control

Excellent long-term power and pulse-to-pulse stability

Three-Channel Wavelength-Tunable Femtosecond Laser




Watt-level output at high repetition rate for fast imaging

Two tunable and one fixed output for simultaneous multibeam excitation


Automated GDD control for shortest pulses at the sample

Industrial-grade design for high power and beam stability


CRONUS-2P tuning curve




CRONUS-2P GDD control range



CRONUS-2P relative intensity noise (RIN)



CRONUS-2P typical output power stability at 900 nm



Specifications

Model	CRONUS-2P		
	Output A	Output B	Output C
Tuning range ¹⁾	680 – 960 nm	940 – 1300 nm	1025 ± 10 nm (fixed)
Output power ²⁾³⁾	> 3 W @ 920 nm	> 2.5 W @ 1100 nm	> 2.5 W
Pulse duration ⁴⁾⁵⁾	< 160 fs		
Repetition rate	77 ± 1 MHz		
Beam quality, M^2 ³⁾⁴⁾	< 1.2		
Polarization	Linear, horizontal		
Beam divergence, full angle	< 1 mrad		< 1.5 mrad
Beam diameter ⁴⁾ ($1/e^2$)	3.0 ± 0.4 mm	3.2 ± 0.4 mm	2.8 ± 0.4 mm
Beam ellipticity ⁴⁾	> 0.8		
Beam astigmatism ⁴⁾	< 20%		
Beam pointing stability ⁶⁾	< 200 µrad		n/a
Long-term power stability, 24 h ⁴⁾⁷⁾	< 1%		
GDD control range	-10 000 to -35 000 fs ² @ 700 nm -3 000 to -20 000 fs ² @ 800 nm 0 to -10 000 fs ² @ 920 nm	0 to -10 000 fs ² @ 960 nm -3 000 to -10 000 fs ² @ 1100 nm -8 000 to -12 000 fs ² @ 1300 nm	n/a

OPTIONAL POWER CONTROL

Tuning range ⁸⁾	680 – 960 nm	940 – 1300 nm	1025 ± 10 nm (fixed)
Output power ⁹⁾	> 2 W @ 920 nm	> 2 W @ 1100 nm	> 1.5 W
Rise/fall time ¹⁰⁾	< 300 ns		
Contrast ratio	1000:1		
GDD control range	0 to -6 500 fs ² @ 920 nm	0 to -10 000 fs ² @ 1100 nm	n/a

OPTIONAL WAVELENGTH EXTENSIONS (UV – VIS)

Second harmonic tuning range	340 – 480 nm ¹¹⁾	480 – 650 nm ¹¹⁾	n/a
Conversion efficiency at peak	> 30%		

ENVIRONMENTAL REQUIREMENTS & DIMENSIONS

Refer to www.lightcon.com

¹⁾ Configuration with dual-output A or dual-output B is also available.

²⁾ Simultaneous mode: > 1 W @ 920 nm, > 1 W @ 1100 nm, and > 2.5 W @ 1025 nm.

³⁾ Power control using AOM is applicable, specifications below.

⁴⁾ Specified at 920 nm, 1100 nm, and 1025 nm, respectively.

⁵⁾ IR pulse duration determined assuming sech2 shape.

⁶⁾ Beam pointing deviation over the entire tuning and GDD control range.


⁷⁾ Expressed as normalized root mean squared deviation (NRMSD); with less than ±1 °C temperature change after 1 h warm up.

⁸⁾ Configuration with dual-output A or dual-output B is also available.

⁹⁾ Simultaneous mode: > 0.7 W @ 920 nm, > 0.7 W @ 1100 nm, and > 1.5 W @ 1025 nm.


¹⁰⁾ Specified from 5% to 95%.

¹¹⁾ Multiple second harmonic configurations available. For more information contact sales@lightcon.com.




Drawings

CRONUS-2P drawing



Laser Source for Advanced Nonlinear Microscopy



High pulse energy for deep imaging


1250 – 1800 nm tuning range for 3P imaging

Down to 50 fs pulse duration for high peak power


Automated wavelength and GDD control for optimal signal

Market-leading pulse-to-pulse energy stability


CRONUS-3P output power and pulse energy vs wavelength, at 1 MHz




CRONUS-3P GDD control range



CRONUS-3P typical pulse duration at 1300 nm and 1700 nm



CRONUS-3P beam profile at 1300 nm



Specifications

Model	CRONUS-3P		CRONUS-3P with power control	
Tuning range	1250 – 1800 nm			
Repetition rate ¹⁾	Single-shot – 1 MHz or 2 MHz			
	1300 nm	1700 nm	1300 nm	1700 nm
Pulse duration	< 50 fs	< 65 fs	< 50 fs	< 65 fs
Output power	> 1100 mW @ 1 MHz > 800 mW @ 2 MHz	> 800 mW @ 1 MHz > 500 mW @ 2 MHz	> 1000 mW @ 1 MHz > 700 mW @ 2 MHz	> 700 mW @ 1 MHz > 400 mW @ 2 MHz
GDD control range ²⁾	-4 000 to +9 000 fs ²	-12 000 to +3 500 fs ²	-4 000 to +9 000 fs ²	-12 000 to +3 500 fs ²
Beam diameter ³⁾	2 – 4 mm			
Beam quality, M ²	< 1.2			
Beam ellipticity	> 0.8			
Beam divergence	< 1 mrad			
Beam pointing stability	< 100 µrad			
Long-term power stability, 24 h ⁴⁾	< 1%			
Pulse-to-pulse energy stability, 1 min ⁴⁾	< 1%			

MAIN OUTPUT WITHOUT GDD CONTROL

Output power ⁵⁾	> 1500 mW @ 1 MHz > 1000 mW @ 2 MHz	> 1050 mW @ 1 MHz > 700 mW @ 2 MHz	n/a
----------------------------	--	---------------------------------------	-----

ADDITIONAL OUTPUTS


Auxiliary 1030 nm amplifier output	1030 ± 10 nm, up to 40 W, up to 2 MHz, < 250 fs
Optional 680 – 920 nm amplifier output	680 – 920 nm, > 1500 mW @ 1 MHz or > 1000 mW @ 2 MHz (@ 920 nm), < 290 fs (compressible to < 50 fs) ⁶⁾
Optional 1030 nm oscillator output	1030 ± 10 nm, up to 500 mW, ≈ 65 MHz, ≈ 200 fs

ENVIRONMENTAL REQUIREMENTS & DIMENSIONS


Refer to www.lightcon.com

- ¹⁾ Lower repetition rate with a higher pulse energy option available.
²⁾ Continuous dispersion control; -4000 fs² compensates a microscope with +4000 fs².
³⁾ 1/e², measured at compressor output.
⁴⁾ Expressed as normalized root mean squared deviation (NRMSD).


- ⁵⁾ Available only for v1. Contact sales@lightcon.com for more details.
⁶⁾ With external compressor without GDD control, < 70% transmission at 920 nm.



CRONUS-3P typical long-term power stability at 1300 nm





CRONUS-3P typical pulse-to-pulse energy distribution at 1300 nm



Drawings

CRONUS-3P drawing





Wavelength-Tunable Sources

LIGHT CONVERSION's OPAs offer a broad tuning range from deep-UV to mid-IR. Coupled with our femtosecond lasers, these OPAs provide an invaluable source for ultrafast spectroscopy, nonlinear microscopy, and a variety of other scientific applications.

I-OPA

The only commercial industrial-grade OPA combines wavelength tunability with robust industrial design.

ORPHEUS | NEO

The next-generation of OPAs with exceptional stability and multiple detectors for continuous power monitoring and diagnostics.

ORPHEUS

Classic OPAs that many are used to. Just like TOPAS, they are quite simple yet offer an extensive range of parameters.

TOPAS

Classic OPAs for Ti:Sapphire lasers.

Continuous wavelength tunability from UV to MIR

Pulse duration from tens of femtoseconds to a few picoseconds

Leading OPA manufacturer for more than 30 years

I-OPA

Industrial-Grade Optical Parametric Amplifier



I-OPA-TW on air-cooled CARBIDE-CB5

Wavelength tunability
in an industrial design


Single-box solution

Tunable or fixed-wavelength
models


Plug-and-play installation and
robust performance

The most compact OPA
in the market

I-OPA-HP typical tuning curves
Pump: 40 W, 400 μ J, 100 kHz



I-OPA-F typical tuning curves
Pump: 40 W, 400 μ J, 100 kHz



Specifications

Model	I-OPA-HP	I-OPA-F	I-OPA-ONE
Configuration	ORPHEUS	ORPHEUS-F	ORPHEUS-ONE
Pump power	Up to 40 W		
Pump pulse energy	20 – 400 μ J		
Repetition rate	Up to 2 MHz		
Tuning range ¹⁾	640 – 1010 nm (signal) 1050 – 2600 nm (idler)	650 – 920 nm (signal) 1200 – 2500 nm (idler)	1350 – 2000 nm (signal) 2100 – 4500 nm (idler)
Conversion efficiency	> 7% @ 700 nm (40 – 400 μ J pump; up to 1 MHz)		> 9% @ 1550 nm (40 – 400 μ J pump; up to 1 MHz)
	> 3.5% @ 700 nm (20 – 40 μ J pump; up to 2 MHz)		> 6% @ 1550 nm (20 – 40 μ J pump; up to 2 MHz)
Spectral bandwidth ²⁾	80 – 220 cm^{-1} @ 700 – 960 nm	200 – 1000 cm^{-1} @ 650 – 920 nm 150 – 1000 cm^{-1} @ 1200 – 2000 nm	60 – 150 cm^{-1} @ 1450 – 2000 nm
Pulse duration ²⁾³⁾	120 – 250 fs	< 55 fs @ 800 – 920 nm < 70 fs @ 650 – 800 nm < 100 fs @ 1200 – 2000 nm	100 – 300 fs
Long-term power stability, 8 h ⁴⁾	< 1% @ 800 nm		< 1% @ 1550 nm
Pulse-to-pulse energy stability, 1 min ⁴⁾	< 1% @ 800 nm		< 1% @ 1550 nm
Wavelength extension options	320 – 505 nm (SHS) ⁵⁾ 525 – 640 nm (SHI) ⁵⁾	Contact sales@lightcon.com	4500 – 10000 nm (DFG)
Pulse compression options ²⁾	n/a	SCMP (signal pulse compressor) ICMP (idler pulse compressor) GDD-CMP (compressor with GDD control)	n/a

PUMP LASER REQUIREMENTS

Pump laser	PHAROS or CARBIDE
Center wavelength	1030 \pm 10 nm
Maximum pump power	40 W
Maximum repetition rate	Up to 2 MHz
Pump pulse energy	20 – 400 μ J
Pulse duration	180 – 300 fs

ENVIRONMENTAL & UTILITY REQUIREMENTS

Operating temperature ⁶⁾	19 – 25 $^{\circ}$ C (air conditioning recommended)
Relative humidity ⁶⁾	20 – 70% (non-condensing)
Electrical requirements	n/a ⁷⁾

¹⁾ In case of fixed wavelength (FW), a single wavelength can be selected from the signal or idler range. The signal may have an accessible idler pair, and vice versa.

²⁾ I-OPA-F broad-bandwidth pulses are compressed externally. Typical pulse duration before compression: 120 – 250 fs, after compression: 25 – 70 fs @ 650 – 900 nm, 40 – 100 fs @ 1200 – 2000 nm.


³⁾ Output pulse duration depends on selected wavelength and pump laser pulse duration.

⁴⁾ Expressed as normalized root mean squared deviation (NRMSD).

⁵⁾ Conversion efficiency is 1.2% at peak; specified as a percentage of pump power.

⁶⁾ Specifications are guaranteed for a maximum temperature variation of \pm 1 $^{\circ}$ C and humidity variation of \pm 10%.

⁷⁾ I-OPA is powered by the same electrical source as the pump laser. Thus, refer to the pump laser electrical requirements.



Broad-Bandwidth Hybrid Optical Parametric Amplifier



Combination of best collinear and non-collinear OPA features

Ultrashort pulses in NIR (650 – 900 nm and 1200 – 2500 nm)


Single-shot – 2 MHz repetition rate

< 100 fs pulse duration


Adjustable spectral bandwidth

Optional long pulse mode for gap-free tunability


ORPHEUS-F typical tuning curves.
Pump: 40 W, 40 μ J, 1000 kHz



ORPHEUS-F typical spectral bandwidth



ORPHEUS-F pulse duration after compression



For custom tuning curves visit
<http://toolbox.lightcon.com/tools/tuningcurves>

Specifications

MAIN OUTPUT (650 – 900 nm and 1200 – 2500 nm)

Mode of operation	Short pulse mode ¹⁾	Long pulse mode
Tuning range	650 – 900 nm (signal) 1200 – 2500 nm (idler)	650 – 1010 nm (signal) 1050 – 2500 nm (idler)
Maximum pump power	80 W	
Pump pulse energy	10 – 500 μJ	
Conversion efficiency ²⁾	> 7% @ 700 nm	
Integrated 2H (515 nm) generation efficiency ³⁾	> 35%	
Pulse duration before compression ¹⁾	< 290 fs	
Spectral bandwidth	200 – 750 cm ⁻¹ @ 650 – 900 nm	60 – 220 cm ⁻¹ @ 650 – 900 nm
Pulse duration after compressor ¹⁾	< 55 fs @ 800 – 900 nm < 70 fs @ 650 – 800 nm < 100 fs @ 1200 – 2000 nm	n/a
Compressor transmission	> 65% @ 650 – 900 nm > 80% @ 1200 – 2000 nm	
Long-term power stability, 8 h ⁴⁾	< 2% @ 800 nm	
Pulse-to-pulse energy stability, 1 min ⁴⁾	< 2% @ 800 nm	

WAVELENGTH EXTENSION OPTIONS (325 – 15000 nm) ⁵⁾

325 – 450 nm (SHS)	> 1%	n/a
325 – 505 nm (SHS)	n/a	> 1%
525 – 650 nm (SHI)		> 0.5%
600 – 650 nm (SHI)	> 0.5%	n/a
210 – 252 nm (FHS)	n/a	> 0.1%
263 – 325 nm (FHI)		
2500 – 15000 nm	See ORPHEUS-MIR;	

PUMP LASER REQUIREMENTS

Pump laser	PHAROS or CARBIDE
Center wavelength	1030 ± 10 nm
Maximum pump power	80 W
Repetition rate	Single-shot – 2 MHz
Pump pulse energy	10 – 500 μJ
Pulse duration ⁶⁾	180 – 500 fs

ENVIRONMENTAL & UTILITY REQUIREMENTS

Refer to www.lightcon.com

¹⁾ In short pulse mode, broadband pulses are compressed externally. Typical pulse duration before compression: 120 – 250 fs, after compression: 25 – 70 fs @ 650 – 900 nm, 40 – 100 fs @ 1200 – 2000 nm.


²⁾ Specified as a percentage of pump power, before compressor. Conversion efficiency at peak is equal to 10% for signal and idler combined.

³⁾ At designated output port; not simultaneous to OPA output.

⁴⁾ Expressed as normalized root mean squared deviation (NRMSD).


⁵⁾ For > 15 μJ pump pulse energy.

⁶⁾ FWHM, assuming Gaussian pulse shape.



Drawings

ORPHEUS-F drawings



Dual Optical Parametric Amplifier



Two simultaneous independently tunable outputs

210 – 16000 nm tuning range

Single-shot – 2 MHz repetition rate

Up to 60 W, 0.5 mJ pump

Compact and cost-effective

CEP-stable option

Specifications

MAIN OUTPUT

Tuning range	Choice between ORPHEUS, ORPHEUS-F, and ORPHEUS-ONE configurations
Output pulse energy	Depends on the configuration, see the specifications of the chosen models
Spectral bandwidth	Depends on configuration, 100 – 750 cm^{-1}
Pulse duration	Depends on configuration, down to < 50 fs
Repetition rate	Single-shot – 2 MHz

PUMP LASER REQUIREMENTS


Pump laser	PHAROS or CARBIDE
Center wavelength	1030 \pm 10 nm
Maximum pump power	60 W
Repetition rate	Single-shot – 2 MHz
Pump pulse energy	16 – 500 μJ
Pulse duration ¹⁾	180-300 fs

ENVIRONMENTAL & UTILITY REQUIREMENTS


Operating temperature ²⁾	19 – 25 °C (air conditioning recommended)
Relative humidity ²⁾	20 – 70% (non-condensing)
Electrical requirements	100 – 240 V AC, 4.5 A; 50 – 60 Hz
Rated power	280 W
Power consumption	Standby: 20 W Max during wavelength tuning: 200 W

¹⁾ FWHM, assuming Gaussian pulse shape.


²⁾ Specifications are guaranteed for a maximum temperature variation of $\pm 1^\circ\text{C}$ and humidity variation of $\pm 10\%$.



ORPHEUS-TWINS (-ONE/-F configuration)
tuning curves.
Pump: 40 W, 40 μ J, 1000 kHz




ORPHEUS-TWINS (ORPHEUS / ORPHEUS configuration)
tuning curves.
Pump: 20 W, 20 μ J, 100 kHz



For custom tuning curves visit
<http://toolbox.lightcon.com/tools/tuningcurves>

Drawings

ORPHEUS-TWINS drawings





Femtosecond Lasers

LIGHT CONVERSION is world-renowned for its industrial-grade Yb-based femtosecond lasers, covering a wide range of scientific, industrial, and medical applications.

FLINT

Expanding the parameter range with repetition rates ranging from 10 to 100 MHz, with power up to 20 W and pulse duration down to 50 fs.

PHAROS

Scientific flexibility and process-tailored output parameters, providing pulse duration down to 100 fs and pulse energy of up to 4 mJ.

CARBIDE

Compact industrial design in air-cooled and water-cooled models, providing up to 120 W, 1 mJ or 80 W, 2 mJ with excellent output stability.

High average power and high pulse energy at a high repetition rate

Market-proven industrial-grade stability and reliability

Tailored to the needs of industry and science

High-Repetition-Rate Lasers



FLINT-FL1

From 10 to 100 MHz repetition rate

Down to 50 fs pulse duration


High-power models, up to 20 W

High-energy energy models, up to 0.5 μ J


Industrial-grade design for high output stability

CEP stabilization or repetition rate locking


FLINT-FL1
Typical spectrum




FLINT-FL1
Typical pulse duration




FLINT-FL1
Typical beam profile




FLINT-FL2-SP
Typical spectrum



FLINT-FL2-SP
Typical pulse duration



FLINT-FL2-SP
Typical beam profile



Specifications

Model	FL1			FL2-SP	FL2		
Key feature	CEP	RRL	Compact	Short pulse	High power and high energy		
Pulse duration	< 100 fs		< 120 fs	< 50 fs	< 120 fs	< 170 fs ¹⁾	
Repetition rate	60 – 100 MHz ²⁾			10 MHz	10 MHz	40 MHz	80 MHz
Maximum output power	0.5 W	1 W	8 W	4 W	5 W	20 W	
Maximum pulse energy	6 nJ ³⁾	12.5 nJ ³⁾	100 nJ ³⁾	0.4 μJ	0.5 μJ		0.25 μJ
Center wavelength	1035 ± 10 nm			1030 ± 10 nm	1030 ± 10 nm		
Polarization	Linear, horizontal						
Beam quality, M ²	< 1.2			< 1.3	< 1.2		
Beam pointing stability	< 10 μrad/°C						
Long-term power stability, 100 h ⁴⁾	< 0.5%						
Integrated 2H generator ⁵⁾	n/a				Optional; conversion efficiency > 30% ⁶⁾		
External 2H, 3H, or 4H generator ⁵⁾	Optional; Refer to HIRO for FLINT						
Integrated attenuator	n/a			Included			

PHYSICAL DIMENSIONS

Laser head (L × W × H)	448 × 206 × 115 mm	543 × 322 × 146 mm
Power supply and chiller rack (L × W × H)	642 × 553 × 540 mm	642 × 553 × 673 mm
Chiller	Different options available. Contact sales@lightcon.com	

ENVIRONMENTAL AND UTILITY REQUIREMENTS

Operating temperature	15 – 30 °C (air conditioning recommended)	
Relative humidity	< 80% (non-condensing)	
Electrical requirements	100 V AC, 7 A – 240 V AC, 3 A; 50 – 60 Hz	100 V AC, 12 A – 240 V AC, 5 A; 50 – 60 Hz
Rated power	200 W	
Power consumption	Laser	100 W
	Chiller	600 W
		150 W
		1000 W

¹⁾ For 20 W output power. Lower power models: 8 W and 12 W, are available upon request.


²⁾ Standard repetition rate is 80 MHz; custom repetition rate can be factory preset from the given range.

³⁾ Depends on the repetition rate. Values are given for 80 MHz.

⁴⁾ With enabled power-lock, under stable environmental conditions. Expressed as normalized root mean squared deviation (NRMSD).


⁵⁾ For external 2H, or even 3H and 4H generation, refer to HIRO for FLINT.

⁶⁾ Conversion efficiency specified at maximum power.




Stability

FLINT-FL2 (20 W) output power stability under harsh environmental conditions over 7 days



FLINT oscillator relative intensity noise (RIN), shot-noise limited at -160 dBc/Hz above 1 MHz



PHAROS

Modular-Design Femtosecond Lasers for Industry and Science



Tunable pulse duration, 100 fs – 20 ps

Maximum pulse energy of up to 4 mJ

Down to < 100 fs right at the output


Pulse-on-demand and BiBurst for pulse control

Up to 5th harmonic or tunable extensions


CEP stabilization or repetition rate locking

Thermally-stabilized and sealed design


PHAROS-PH2-UP
Typical pulse duration



PHAROS-PH2-UP
Typical spectrum



PHAROS
Pulse energy vs fundamental repetition rate



Specifications

Model	PH2-10W	PH2-20W-SP			PH2-4mJ	PH2-UP	
-------	---------	------------	--	--	---------	--------	--

OUTPUT CHARACTERISTICS

Center wavelength ¹⁾	1030 ± 10 nm						
Maximum output power	10 W	20 W					
Pulse duration ²⁾	< 290 fs	< 190 fs			< 450 fs ³⁾	< 100 fs	
Pulse duration tuning range	290 fs – 10 ps (20 ps on request)	190 fs – 10 ps (20 ps on request)			450 fs – 10 ps	100 fs – 10 ps	
Maximum pulse energy	0.2 mJ	0.4 mJ	1 mJ	2 mJ	4 mJ	0.4 mJ	1 mJ
Repetition rate	Single-shot – 1 MHz						
Pulse selection	Single-shot, pulse-on-demand, any fundamental repetition rate division						
Polarization	Linear, horizontal						
Beam quality, M ²	< 1.2	< 1.3				< 1.2	
Beam diameter ⁴⁾	3.3 ± 0.5 mm	4.0 ± 0.5 mm	4.5 ± 0.5 mm	6.8 ± 0.7 mm		4.5 ± 0.5 mm	6 ± 0.5 mm
Beam pointing stability	< 20 µrad/°C						
Pre-pulse contrast	< 1:1000						
Post-pulse contrast	< 1:200						
Pulse-to-pulse energy stability, 24 h ⁵⁾	< 0.5%						
Long-term power stability, 100 h ⁵⁾	< 0.5%						

MAIN OPTIONS

Oscillator output ⁶⁾	1 – 7 W, 50 – 250 fs, ≈ 1035 nm, ≈ 76 MHz						
Harmonic generator ⁷⁾	515 nm, 343 nm, 257 nm, or 206 nm						
Optical parametric amplifier ⁸⁾	320 – 10000 nm; see page 12						
BiBurst option	Tunable GHz and MHz burst with burst-in-burst capability						
CEP stabilization	Optional						
Repetition rate locking							

PHYSICAL DIMENSIONS

Laser head (L × W × H) ⁹⁾	730 × 419 × 230 mm	827 × 492 × 250 mm	770 × 419 × 230 mm
Chiller (L × W × H)	590 × 484 × 267 mm		
24 V DC power supply (L × W × H) ⁹⁾	280 × 144 × 49 mm		

ENVIRONMENTAL & UTILITY REQUIREMENTS

Operating temperature	15 – 30 °C (air conditioning recommended)		
Relative humidity	< 80% (non-condensing)		
Electrical requirements	Laser	100 V AC, 12 A – 240 V AC, 5 A, 50 – 60 Hz	
	Chiller	100 – 230 V AC, 50 – 60 Hz	
Rated power	Laser	1000 W	
	Chiller	1400 W	
Power consumption	Laser	600 W	
	Chiller	1000 W	

¹⁾ Precise wavelengths for specific models are available upon request.

²⁾ Assuming Gaussian pulse shape.

³⁾ Pulse duration can be reduced to < 250 fs if pulse peak intensity of > 50 GW/cm² is tolerated by the customer setup.

⁴⁾ FW 1/e², measured at laser output, using maximum pulse energy.


⁵⁾ Under stable environmental conditions. Expressed as normalized root mean squared deviation (NRMDS).

⁶⁾ Available simultaneously. Contact sales@lightcon.com for more details or customized solutions.

⁷⁾ Integrated. For external harmonic generator, refer to HIRO.

⁸⁾ Integrated. For more options and OPAs for -4mJ and -UP models, refer to www.lightcon.com.


⁹⁾ Dimensions depend on laser configuration and integrated options.



Beam properties


PHAROS

Typical M^2 measurement data




PHAROS

Typical near-field beam profile



PHAROS


Typical far-field beam profile




Stability measurements

PHAROS


Long-term power stability



Output power of industrial-grade PHAROS lasers operating 24/7 and the current of the pump diodes over the years




PHAROS output power and the stability of beam direction with power lock enabled, across varying environmental conditions



PHAROS

Typical pulse-to-pulse energy stability




CEP stabilization

PHAROS lasers can be equipped with feedback electronics for carrier-envelope phase (CEP) stabilization of the output pulses. The carrier-envelope offset (CEO) of the PHAROS oscillator is actively locked to 1/4th of the repetition rate with a < 100 mrad standard deviation. The CEP stable pulses from the

synchronized amplifier have a < 350 mrad standard deviation. The CEP drift occurring inside the amplifier and the user's setup can be compensated with an out of loop f-2f interferometer, which is a part of the complete PHAROS active CEP stabilization package.


PHAROS

Short-term CEP stability operating at 200 kHz repetition rate



PHAROS

Long-term CEP stability operating at 200 kHz repetition rate




Repetition rate locking

The oscillator of PHAROS laser can be customized for repetition rate locking applications. Coupled with the necessary feedback electronics, the repetition rate is synchronized to an external RF source using the two piezo stages installed inside the cavity.


The repetition rate locking system can assure an integrated timing jitter of less than 200 fs for RF reference frequencies larger than 500 MHz. Continuous phase shifting is available on request.

Phase noise data of PHAROS oscillator locked to a 2.8 GHz RF source



Timing jitter stability over 14 h


PHAROS oscillator locked to a 2.8 GHz RF source



Drawings


PHAROS-PH2-730 drawing.

PH2 or PH2-SP with FEC, BiBurst, or harmonics




PHAROS-PH2-827 drawing

PH2 with -HE harmonics, PH2-4mJ, or PH2-UP with harmonics



Unibody-Design Femtosecond Lasers for Industry and Science



CARBIDE-CB3

Tunable pulse duration,
190 fs – 20 ps

Maximum output of
120 W, 1 mJ or 80 W, 2 mJ

Single-shot – 10 MHz
repetition rate

NEW


Pulse-on-demand and
BiBurst for pulse control

Up to 5th harmonic or
tunable extensions


Air-cooled and
water-cooled models

Compact industrial-grade design


CARBIDE-CB3
Typical pulse duration




CARBIDE-CB3
Typical spectrum



CARBIDE-CB3
Typical beam profile



CARBIDE-CB3-120W
Long-term power stability



CARBIDE-CB3 specifications

NEW

Model	CB3-20W	CB3-40W	CB3-80W	CB3-120W
-------	---------	---------	---------	----------

OUTPUT CHARACTERISTICS

Cooling method	Water-cooled			
Center wavelength	1030 ± 10 nm			
Maximum output power	20 W	40 W	80 W	120 W
Pulse duration ¹⁾	< 250 fs		< 350 fs ²⁾	< 250 fs
Pulse duration tuning range	250 fs – 10 ps		350 fs – 10 ps	250 fs – 10 ps
Maximum pulse energy	0.4 mJ	0.2 mJ	0.8 mJ	2 mJ
Repetition rate	Single-shot – 1 MHz	Single-shot-1 MHz (2 MHz on request)	Single-shot – 10 MHz	Single-shot – 2 MHz
Pulse selection	Single-shot, pulse-on-demand, any fundamental repetition rate division			
Polarization	Linear, vertical; 1: 1000			
Beam quality, M ²	< 1.2			
Beam diameter ³⁾	3.9 ± 0.4 mm		4.2 ± 0.4 mm	5.1 ± 0.7 mm
Beam pointing stability	< 20 μrad/°C			
Pulse energy control	FEC ⁴⁾	Attenuator ⁵⁾	FEC ⁴⁾	
Pulse picker leakage	< 0.25%	< 0.5%	< 0.25%	
Pulse-to-pulse energy stability, 24 h ⁶⁾	< 0.5%			
Long-term power stability, 100 h ⁶⁾	< 0.5%			

MAIN OPTIONS

Oscillator output	< 0.5 W, 120 – 250 fs, 1030 ± 10 nm, ≈ 65 MHz ⁷⁾			
Harmonic generator ⁸⁾	515 nm, 343 nm, 257 nm, or 206 nm			
Optical parametric amplifier ⁹⁾	320 – 10000 nm; see page 12			n/a
BiBurst option	Tunable GHz and MHz burst with burst-in-burst capability			

PHYSICAL DIMENSIONS

Laser head (L × W × H)	633 × 350 × 174 mm		
Chiller (L × W × H)	585 × 484 × 221 mm	680 × 484 × 307 mm	
24 V DC power supply (L × W × H)	280 × 144 × 49 mm ¹⁰⁾	320 × 200 × 75 mm	376 × 449 × 88 mm

ENVIRONMENTAL AND UTILITY REQUIREMENTS

Operating temperature	15 – 30 °C		
Relative humidity	< 80% (non-condensing)		
Electrical requirements	Laser	100 V AC, 7 A – 240 V AC, 3A; 50 – 60 Hz	100 V AC, 12 A – 240 V AC, 5 A 50 – 60 Hz
	Chiller	100 – 230 V AC; 50 – 60 Hz	200 – 230 V AC; 50 – 60 Hz
Rated power	Laser	600 W	1000 W
	Chiller	1400 W	2000 W
Power consumption	Laser	500 W	900 W
	Chiller	1000 W	1300 W

¹⁾ Assuming Gaussian pulse shape.

²⁾ Pulse duration can be reduced to < 250 fs if pulse peak intensity of > 50 GW/cm² is tolerated by the customer setup.

³⁾ FW 1/e², using maximum pulse energy.

⁴⁾ Fast energy control (FEC) provides fast, full-scale individual pulse energy control; an external analog control input is available.

⁵⁾ Waveplate-based variable optical attenuator (VOA); an external analog control input is available.


⁶⁾ Under stable environmental conditions. Expressed as normalized root mean squared deviation (NRMSD)

⁷⁾ Available simultaneously, requires a scientific interface. Contact sales@lightcon.com for more details or customized solutions.

⁸⁾ Integrated. For external harmonic generator, refer to HIRO.

⁹⁾ Integrated. For more options and OPAs, refer to www.lightcon.com.

¹⁰⁾ Power supply can be different if optional 2 MHz version is selected.



CARBIDE-CB5 (air-cooled) specifications

Model	CB5		CB5-SP
OUTPUT CHARACTERISTICS			
Cooling method	Air-cooled ¹⁾		
Center wavelength	1030 ± 10 nm		
Maximum output power	6 W	5 W	
Pulse duration ²⁾	< 290 fs		< 190 fs
Pulse duration tuning range	290 fs – 20 ps		190 fs – 20 ps
Maximum pulse energy	100 µJ	83 µJ	100 µJ
Repetition rate	Single-shot – 1 MHz		
Pulse selection	Single-shot, pulse-on-demand, any fundamental repetition rate division		
Polarization	Linear, vertical; 1: 1000		
Beam quality, M ²	< 1.2		
Beam diameter ³⁾	2.1 ± 0.4 mm		
Beam pointing stability	< 20 µrad/°C		
Pulse energy control	Attenuator ⁴⁾	AOM ⁵⁾	Attenuator ⁴⁾
Pulse picker leakage	< 2%	< 0.1%	< 2%
Pulse-to-pulse energy stability, 24 h ⁶⁾	< 0.5%		
Long-term power stability, 100 h ⁶⁾	< 0.5%		
MAIN OPTIONS			
Oscillator output	n/a		
Harmonic generator ⁷⁾	515 nm, 343 nm, 257 nm, or 206 nm;		
Optical parametric amplifier ⁸⁾	320 – 10000 nm; see page 12		
BiBurst option	n/a		
PHYSICAL DIMENSIONS			
Laser head (L × W × H)	633 × 324 × 162 mm		
Chiller	Not required		
24 V DC power supply (L × W × H)	220 × 95 × 46 mm		
ENVIRONMENTAL AND UTILITY REQUIREMENTS			
Operating temperature	17 – 27 °C		
Relative humidity	< 80% (non-condensing)		
Electrical requirements	100 V AC, 3 A – 240 V AC, 1.3 A; 50 – 60 Hz		
Rated power	300 W		
Power consumption	150 W		

¹⁾ Water-cooled version available on request.

²⁾ Assuming Gaussian pulse shape.

³⁾ $FW\ 1/e^2$, using maximum pulse energy.


⁴⁾ Waveplate-based variable optical attenuator (VOA); an external analog control input is available.

⁵⁾ Enhanced contrast AOM. Provides fast amplitude control of output pulse train.


⁶⁾ Under stable environmental conditions. Expressed as normalized root mean squared deviation (NRMSD).

⁷⁾ Integrated. For external harmonic generator, refer to HIRO.


⁸⁾ Integrated. For stand-alone OPAs, refer to www.lightcon.com.




CARBIDE-CB5-SP
Typical pulse duration



CARBIDE-CB5
Long-term power stability




CARBIDE-CB5
Typical beam profile




Stability measurements

CARBIDE-CB3 output power and beam direction stability with power lock enabled, across varying environmental conditions




CARBIDE-CB3
Typical pulse-to-pulse energy stability





Drawings

CARBIDE-CB3 drawing



Air-cooled CARBIDE-CB5 with attenuator drawing





Harmonic imaging of the heart's conductive system using FLINT femtosecond oscillator.
Courtesy of the Biomedical Photonics Laboratory, Vilnius University.

Nonlinear Microscopy Applications

LIGHT CONVERSION delivers best-in-class lasers and laser systems for today's most demanding applications.

Functional 3P neuroimaging



2P optogenetics

SHG, THG, and 2P imaging

Functional 3P neuroimaging

Recording real-time single-neuron activity in the deep brain layers of awake animals is essential for understanding behavior, brain connectivity, and function. These applications have been advanced by neuron imaging and stimulation techniques using high-power, high-pulse-energy lasers with medium-repetition rates, tunable in the SWIR range, which aligns with the biological transparency

windows at 1300 nm and 1700 nm. For 2P and 3P excited fluorescence, and harmonic-generation (SHG, THG) imaging in deep tissues, dispersion-controlled femtosecond pulses from I-OPA and ORPHEUS OPAs and microscopy-dedicated CRONUS lasers represent state-of-the-art choices.




3P microscopy with adaptive optics for focus sensing and shaping to compensate for both aberrations and scattering. ORPHEUS-F excitation at 1300 nm enabled imaging up to 1.1 mm below the pia within the intact brain.

Courtesy of Jianan Y. Qu group, the Hong Kong University of Science and Technology. Source: Qin et al., Deep tissue multi-photon imaging using adaptive optics with direct focus sensing and shaping, Nature Biotechnology 40 (2022).

2P and 3P calcium imaging at depth in mouse brain

Three-photon microscopy (3PM) has gained popularity as a tool able to extend the capabilities of two-photon microscopy (2PM) by imaging deeper layers in the brain and other tissues such as tumors and bone.

Imaging depth in 2PM is limited by the scattering and absorption of excitation light within the tissue. 3PM overcomes this limit because the higher nonlinearity of the 3P excitation reduces the background.






Comparison of in vivo 2P and 3P calcium imaging of mouse visual cortex GCaMP neurons on a Thorlabs Bergamo II microscope using a typical 2P laser and Light Conversion's CRONUS-3P (3P) laser at 920 nm and 1300 nm, respectively.

Courtesy of CSHL ISFNS 2024 school organizers, Willis Broden Jr. and Sergey Matveev (Thorlabs).

2P optogenetics

Despite the advances in 3-photon excitation sources providing longer wavelengths and higher pulse energies, certain imaging challenges are still better addressed by tunable high-repetition-rate oscillator-based lasers. This is especially true when imaging speed is the primary factor. For these applications, the **CRONUS-2P** laser offers the ultimate solution with its optically synchronized three

outputs, two of which are independently tunable. A three-beam source enables a variety of multiphoton excitation pathways, many of which are inaccessible using traditional single- and two-beam solutions. Furthermore, the independent tunability of the two beams enables new coherent Raman scattering modalities.




2P optogenetic stimulation of individual neurons using **CRONUS-2P**.

Courtesy of Albrecht Stroth group, University Medical Center Mainz and Leibniz Institute for Resilience Research. Source: T. Fu et al., Exploring two-photon optogenetics beyond 1100 nm for specific and effective all-optical physiology, *iScience* 24 (2021).

Raster-scanning 2P/3P microscopy

For applications requiring a fixed-wavelength femtosecond laser, such as multiphoton-driven fluorescence, excited at 1 μm , and harmonic-generation (SHG, THG) microscopy, the **FLINT** oscillator is a high-performance solid-state source in a proven, industrial-

grade package and a compact footprint. The **FLINT** oscillator provides stable 24/7 operation with excellent noise performance, characterized by a RIN of < 140 dBc/Hz above 200 kHz and shot-noise-limited performance at -160 dBc/Hz above 1 MHz.



SHG and THG images of H&E-stained colon using the **FLINT** femtosecond oscillator.


Courtesy of Virginijus Barzda group, Vilnius University.



Widefield polarimetric SHG microscopy

Cancer diagnosis and surgical treatment rely on imaging techniques that demand specificity and high throughput. Polarization-resolved second-harmonic generation (P-SHG) microscopy shows potential for visualizing structural changes in collagen networks and the extracellular matrix associated with tumor development. Moreover, P-SHG imaging is label-free and compatible with live tissue imaging at depth. However, traditional raster scanning methods are too slow for clinical applications, and interpreting the structural sensitivity of P-SHG can be challenging.

Nonlinear widefield microscopy addresses these limitations by utilizing amplified femtosecond lasers to increase imaging throughput and field of view. Additionally, machine learning (ML) techniques enable data-driven analysis, facilitating tasks such as automating tumor margin delineation and scoring. By leveraging **PHAROS** and **CARBIDE** lasers in conjunction with ML-augmented widefield microscopy, we can potentially extend the benefits of nonlinear microscopy to the scale required for biomedical and clinical applications.




Large-area widefield P-SHG microscopy of human lung tissue tumor margins conducted using the **PHAROS** laser. Image parameters, including SHG intensity, R-ratio, and degree of circular polarization, as well as SHG circular and linear dichroism, are employed in unsupervised ML algorithms to determine the tumor boundary.

Courtesy of Virginijus Barzda group, University of Toronto, and Brian C. Wilson group, Princess Margaret Cancer Centre. Source: Mirsanaye et al., Unsupervised determination of lung tumor margin with widefield polarimetric second-harmonic generation microscopy, *Scientific Reports* 12 (2022).

SHG, THG, and 2P imaging

Fixation methods, such as formalin, are commonly used for tissue preservation to maintain their structure as close as possible to the native condition. However, these fixatives chemically interact with tissue molecules, potentially altering their structure. To assess the impact of preservation methods, such as chemical fixatives, on

the nonlinear capabilities of protein components within mouse tissues, nonlinear two-photon (2P) microscopy and the **CRONUS-2P** femtosecond laser were utilized. These techniques take advantage of the SHG and THG emission properties of tissue components.




SHG signals from collagen, 2P excitation microscopy and THG signals from elastin in vibratome sections of mouse kidney after different treatments, registered using the **CRONUS-2P** femtosecond laser source.

Courtesy of Frauke Alves and Fernanda Ramos-Gomes, Max-Planck Institute for Multidisciplinary Sciences, Germany.

Combined SHG and THG imaging

Adult zebrafish heart ventricle section used in a scar formation study imaged with the FLINT femtosecond oscillator. The brightfield image is stained with Masson's trichrome (MT), where connective tissue appears blue and muscle appears red/brown.

SHG and THG images reveal collagen and muscle structure at the periphery of the bulbus arteriosus, while MT-stained elastin is visualized in the center in THG.




Adult zebrafish heart ventricle section imaged using the FLINT femtosecond oscillator.

Samples courtesy of Justas Lazutka at the Vilnius University Life Sciences Center. Nonlinear imaging courtesy of the Virginijus Barzda group at the Vilnius University Faculty of Physics.

Label-free in vivo imaging

Understanding biological complexity requires minimally disruptive imaging tools capable of providing multiplexed molecular contrasts. To address this need, S. You's laboratory at the Massachusetts Institute of Technology is developing a non-invasive, label-free microscopy approach using CRONUS-3P to visualize biosystems.

As part of a study on neuropathic pain, the image reveals the rich microenvironment of an unprocessed, intact mouse whisker pad: collagen capsule (green), comprising the follicle with muscles (yellow) supporting it, adipocytes (purple), stromal cells, and immune cells.



Mouse whisker pad using label-free microscopy.

Courtesy of Sixian You group, Massachusetts Institute of Technology.


Global Representative Network

AUSTRALIA NEW ZEALAND	Lastek Pty Ltd. Adelaide, Australia Phone: +61 8 84 438 668 ricardas@lastek.com.au www.lastek.com.au	ISRAEL	ROSH Electroptics Ltd. Natanya, Israel Phone: +972 (0)9 862 7401 info@roshelop.co.il www.roshelop.co.il
BELGIUM, NETHERLANDS, LUXEMBOURG	Laser 2000 Benelux C.V. Vinkeveen, Netherlands Phone: +31 297 266191 info@laser2000.nl www.laser2000.nl	ITALY	Optoprim S.r.l. Vimercate, Italy Phone: +39 039 834 977 info@optoprim.it www.optoprim.it
BRAZIL	Photonics Ltda São Paulo, Brazil Phone: +55 11 2839 3209 info@photonics.com.br www.photonics.com.br	JAPAN	Phototechnica Corp. Saitama, Japan Phone: +81 48 871 0067 voc@phototechnica.co.jp www.phototechnica.co.jp
CZECH REPUBLIC, SLOVAKIA	Femtonika s.r.o. Zbýšov, Czech Republic Phone: +420 792 417 400 info@femtonika.cz www.femtonika.cz	KOREA	Light Conversion Korea Daejeon, Korea Phone: +82 42 368 1010 jungsik.seo@lightcon.com
CHINA	Light Conversion China Shenzhen, China Phone: +86 189 4874 5558 sales.china@cn.lightcon.com	POLAND	Amecam Warszawa, Poland Phone: +48 602 500 680 amecam@amecam.pl www.amecam.pl
	Beijing Light-Quantum Technology Co., Ltd. Beijing, China Phone: +86 10 8290 0415 sales@light-quantum.cn www.light-quantum.cn	SINGAPORE	Acexon Technologies Pte Ltd. Singapore Phone: +65 6565 7300 sales@acexon.com www.acexon.com
	Genuine Optronics Limited Shanghai, China Phone: +86 21 64 325 169 jye@gen-opt.com www.gen-opt.com	SPAIN, PORTUGAL	Innova Scientific S.L. Las Rozas de Madrid, Spain Phone: +34 91 710 56 50 rafael.pereira@innovasci.com www.innovasci.com
FRANCE, SWITZERLAND, BELGIUM	Jean-François Poisson Industrial Market Development Manager Phone: +33 674 48 0778 jf.poisson@lightcon.com	SWITZERLAND	GMP SA Renens, Switzerland Phone: +41 21 633 21 21 info@gmp.ch www.gmp.ch
FRANCE	Frédéric Berthillier Scientific Market Development Manager Phone: +33 745 014 410 frederic.berthillier@lightcon.com	TAIWAN	Alaser Co. Ltd. Taipei, Taiwan Phone: +886 2 2377 3118 alexfu@alaser.com.tw www.alaser.com.tw
GERMANY, AUSTRIA, SWITZERLAND	Ulrich Höchner Industrial Market Development Manager Phone: +49 157 8202 5058 u.hoechner@lightcon.com	TURKEY	Innova Teknoloji Ltd. İstanbul, Turkey Phone: +90 216 315 03 36 eryetistir@innova-teknoloji.com www.innova-teknoloji.com
GERMANY, AUSTRIA	Christian Hellwig Scientific Market Development Manager Phone: +49 174 204 9053 christian.hellwig@lightcon.com	UNITED KINGDOM	Photonic Solutions Ltd. Edinburgh, UK Phone: +44 (0) 131 664 8122 ben.agate@photonicsolutions.co.uk www.photonicsolutions.co.uk
	Stefan Piontek Scientific Market Development Manager Mobile +49 176 8345 7119 stefan.piontek@lightcon.com	USA, CANADA	Light Conversion-USA, Inc. Bozeman, MT, USA Phone: +1 833 685 2872 saleslc@lightcon-usa.com
INDIA	Anatech Laser Instruments Pvt. Ltd. Mumbai, India Phone: +91 22 4121 0001 / 02 / 03 sales@anatechlaser.com www.anatechlaser.com		

Lost in calculations?

Try interactive calculators for scientists & engineers

Optics Toolbox



toolbox.lightcon.com

All calculators on the same page

Optical table layout planner

Custom OPA tuning curves



LIGHT CONVERSION CHINA

702-1, F1 Building, TCL Science Park, No.1001
Zhongshanyuan Road, Nanshan Dist., Shenzhen, China
Phone: +86 189 4874 5558
sales.china@cn.lightcon.com

LIGHT CONVERSION KOREA

520-ho, 65, Techno 3-ro,
Yuseong-gu, Daejeon, 34016, Korea
Phone: +82 42 368 1010
jungsik.seo@lightcon.com

LIGHT CONVERSION

Keramiku 2B, LT-10233 Vilnius, Lithuania

Tel.: +370 5 2491830
Website: www.lightcon.com
Sales: sales@lightcon.com
Service: service@lightcon.com

LIGHT CONVERSION USA, INC

619 N Church Ave, Unit #3
Bozeman, MT 59715
Phone: +1 833 685 2872
Fax: +1 833 395 2872
saleslo@lightcon-usa.com



toolbox.lightcon.com